

Certificate of Analysis

Low Gold Material - *AuBlank1*

Product contains less than 0.002 ppm gold.

The above value applies only to product in sachets that have an identification number of either 78 499 or 84 537.

Finely pulverized material that had previously been analysed and found to have a very low gold content was placed into sachets in a gold-free environment. A total of 1454 sachets were filled. 34 sachets (equivalent to three times the cube root of 1454) were selected at regular intervals from the packaging run and were submitted to a laboratory that specializes in low-level analysis for gold. The samples were analysed by aqua regia digestion followed by enhanced graphite furnace atomic absorption spectroscopy. The average of the 34 analyses was 0.5 parts per billion gold, with the highest value being 1.1 parts per billion gold.

A copy of the analytical report from the independent laboratory is available on request.