



BENCH TOP RING MILL

WHISPER SERIES

USES

Used for pulverising rocks, soil, coal, cement, glass, limestone, bricks, wood, plant material, slags, concrete, etc. for XRF, XRD or other instrumental methods or wet chemistry.

AT A GLANCE

Rocklabs Whisper Series Bench Top Ring Mills have been developed for laboratories pulverising 1 to 100 gram samples.

- Small and compact to fit on a work bench
- Heads/Bowls can also be used on the normal size ring mill
- Noise level less than 80 dBA.
- Motor: 0.18kW

Depth 420 mm
Width 420 mm
Height 500 mm



HEAD OPTIONS AVAILABLE

- | | |
|----------------|--------------------|
| 1 Chrome Steel | 2 Chrome Steel, |
| 3 Carbon steel | with split discus |
| 4 Sialon | 5 Tungsten Carbide |
| 6 Zirconia | |



ADVANTAGES

- Very fast grinding
- Dust free grinding
- Uniformly fine grinding
- Very low contamination of samples
- Easy to clean
- Low capital cost
- Long life
- Maintenance free