



ROCKLABS

WORLD LEADERS IN SAMPLE PREPARATION EQUIPMENT AND REFERENCE MATERIALS FOR USE IN GOLD ASSAYING

March 2008

New Product Release - SMART BOYD-RSD COMBO



The first Smart Boyd-RSD Combo was put into service in 2006 and incorporates all the unique features of the Boyd-RSD Combo. Features include:

- Ø the incoming sample is weighed on a platform balance
- Ø the weight is communicated to the PLC
- Ø the split from the RSD is then adjusted automatically

BOYD-RSD COMBO introduced in 1996



Boyd-RSD incorporates all the unique features of the Boyd Crusher with the addition of a Rotating Sample Divider for “one step splitting and crushing”.

The sample is loaded into the Boyd and in one simple operation, it is crushed, conveyed to the RSD via a vibratory feeder, then split into the proportions of your choice.

Can be upgraded to a Smart Boyd-RSD Combo

BOYD CRUSHER introduced in 1993



Boyd Crusher is a unique patented design of double acting jaws – both top and bottom jaws are driven

Maximum lump size of 55mm can be crushed down to 2mm or less in one pass with feed rate of 5kg per minute.

Can be upgraded to Boyd-RSD Combo or Smart Boyd-RSD Combo.

For further information, please contact your local agent listed on our website or email us.